

# USING EVZIO® AUTO-INJECTOR

- Use EVZIO® exactly as prescribed by your healthcare provider and as described in the *Patient Instructions for Use*.
- Each EVZIO® auto-injector contains only one dose of medicine.
- As part of your opioid emergency plan, practice using the Trainer for EVZIO®.
- Make sure your family and caregiver also practice using the Trainer for EVZIO®.



**STOP**   
**OVERDOSE**

# ADMINISTRATION STEPS

---

- 1 Pull EVZIO<sup>®</sup> from the outer case.
- 2 Pull off the red safety guard.
- 3 Place the black end against the middle of the patient's outer thigh, through clothing (pants, jeans, etc.) if necessary, then press firmly and hold in place for 5 seconds.
- 4 After using EVZIO<sup>®</sup>, get emergency medical help immediately.

For more information, go to  
[evzio.com/patient/about-evzio/index.php](http://evzio.com/patient/about-evzio/index.php)